



CytoCHECK SPAchip[®]

CALCIUM GREEN SINGLE- DETECTION KIT

Instructions for Use

Product Reference: S-002-CAG

For use with flow cytometers and cell imaging platforms

© 2022. **Arrays for Cell Nanodevices SL**

CAIT - UPM. Campus de Montegancedo s/n, 28223 Pozuelo de Alarcón - Madrid (SPAIN)

+34 910 679 519- info@a4cell.com

www.a4cell.com

For research use only

CytoCHECK SPAchip[®] Calcium Single-Detection Kit

Product Reference: S-002-CAG

For use with flow cytometers and cell imaging platforms

1. General Introduction

CytoCHECK SPAchip[®] assay kits are novel cell-based assays for living single-cell developed by A4Cell that bring together the fields of nanotechnology and cell biology. CytoCHECK SPAchip[®] kits are composed of fluorescently labeled silicon microparticles -SPAchips[®]- that can be internalized in the cytosol of cultured cells to monitor changes in specific intracellular analyte concentrations for long periods of time.

CytoCHECK SPAchip[®] Calcium Green Single-Detection Kit allows measurement of intracellular calcium levels by changes in fluorescence intensity, which allows a more comprehensive study of the living single-cell physiology and maximizes the performance of imaging platforms and flow cytometers.

CytoCHECK SPAchip[®] Single-Detection kits are optimized for its use with confocal microscopes and HCS/HCA analyzers with 20X or over magnification objectives, yet epi-fluorescence microscopes and imaging system with fixed wavelength filters can be also used¹. Flow cytometers have additionally been validated to analyze intra and extracellular SPAchips[®].

SPAchip[®] kits can be handled as common fluorophores and chemical probes for cell biology. After resuspending them in the assay buffer (a SPAchip[®]-to-cell ratio of 2:1 is recommended), an overnight incubation is required to allow SPAchip[®] to be incorporated in the cytosol². SPAchip[®] will remain in the cytosol for longer than one month to monitor the evolution of the culture.




! **Storage Information**

Upon receipt, each kit can be stored protected from light at 2-8°C.

¹For using CytoCHECK SPAchip[®] kits in non-confocal imaging systems,

² CytoCHECK SPAchip[®] kits are optimized for its use with adherent cells. Some cells may need a lipofection step to maximize SPAchip internalization. Please contact with A4Cell for further information or questions about your cell model.

CytoCHECK SPChip[®] Calcium Single-Detection Kit Workflow

	<i>SPChip[®] Assay Workflow</i>	<i>Notes</i>	
<i>ASSAY SPChip[®] Preparation</i>	<i>Dissolve ASSAY SPChip[®] tube</i>	Thoroughly dissolve ASSAY SPChip [®] tube solid film in assay buffer (vortexing might be required)	
	<i>Wash with Assay Buffer</i> } x2	Centrifugate tubes at. 4300xg for 10 min. Discard supernatant and keep SPChip [®] containing pellet.	
	<i>Resuspend</i>	Resuspend the SPChips in 1mL ASSAY Buffer Repeat x2. Resuspend the ASSAY SPChips in 100 µL ASSAY Buffer to obtain a stock solution of 2.5x10 ⁴ SPChips/µL	
	<i>Store</i>	ASSAY SPChip [®] stock solution can be stored protected from light at 4°C for up to 6 months.	
<i>Cell Assay</i>	<i>Cell Culture</i>	Seed cells in a multi-well plate. Incubate the cells up to 50-70% confluence ³ .	
	<i>SPChip Dilution</i>	Dilute CONTROL and ASSAY SPChips in cell culture medium to obtain a final SPChip [®] to cell ratio of 2:1.	
	<i>Addition to Cell Culture</i>	Aspire the cell medium and add SPChip [®] -loaded fresh culture medium.	
	<i>Incubation</i>	Incubate overnight in a cell incubator (37°C, 5% CO ₂) to allow SPChip [®] internalization by the cells.	
<i>Fluorescence Quantification</i>	<i>Data Acquisition</i>	Use green fluorescence channel to image SPChip [®] detection (excitation 488 nm, emission 520 nm)	
	<i>Data Analysis</i>		

³ This cell culture protocol only provides a guideline, and it should be modified depending on specific needs. Please, read the whole protocol before starting.

2. Contents and Storage

Each CytoCHECK SPACHIP® Calcium Single-Detection Kit contains:

Materials		Units
ASSAY SPACHIP® tube (embedded in a solid fluorescence-protective soluble film)	~2.5x10 ⁶ ASSAY SPACHIPS	1
ASSAY buffer tube (Sterile, cell culture suitable)	5 mL	1
CONTROL SPACHIP® tube (non-fluorescent, ready-to-use)	~2.5x10 ⁵ CONTROL SPACHIPS/100 uL	1




Figure 1: a) Emission spectra of Calcium--sensing probes when excited at the indicated wavelength. **b)** Brightfield image of SPACHIPS once in solution (Image from Torras et al. 2015 DOI: 10.1002/adma.201504164) Scale bar = 10 μ m.

3. Materials to be Supplied by the User

- Pipettes and pipette tips (1-10 μ L, 2-20 μ L, 20-200 μ L, 200-1000 μ L)
- Multi-channel pipettes and pipette tips (50-300 μ L) (not mandatory but desirable)
- Vortex mixer
- Mini-centrifuge
- Reagent reservoirs
- Cell culture plate (multi-well)
- Cell culture facilities
- Cell culture media (phenol red free is recommended)
- Cell imaging system (HCS/HCS, Confocal microscope, imaging plate reader, etc...). At least 20x magnification objectives are required for quantitative analysis.
- Image analysis software (Contact A4Cell Staff for support)

If quantitative analysis is required, it will be necessary:

- Calcium flux positive and negative controls (for example, A23187 and BAPTA-AM ...)

4. Assay Procedure

NOTE: Work under sterile conditions. Protect the tubes from light, specially from UV light.

a. ASSAY SPChip[®] Preparation

1- Add 1 mL of assay buffer to ASSAY SPChip[®] tubes and mix until complete solubilization of the membrane (vortexing required)

2- Centrifugate the tubes at approx. 4300 xg (8000 RPM in a 6-cm minispin rotor) for 10 minutes.

3-Carefully, aspire and discard supernatant and wash pellets with 1 mL of solubilization medium. Repeat steps 1 and 2.

4-Carefully, aspire and discard supernatant (**Look out not to aspire the pellet!!!**). Resuspend the pellet in **100 µL of assay buffer** to obtain a SPAchip[®] stock solution. This results in approximately **2.5x10⁴ SPAchips/µl**. Once prepared, Assay SPAchip[®] tube can be stored at 2-8°C protected from light for up to 6 months.

b. ASSAY PROTOCOL (standard protocol for 96-well plate with adherent cells)⁴:

NOTE: ASSAY and CONTROL SPAchip[®] dilutions suggested in this section are suggested for 96-well plates. Optimization might be required for different cell types or multiwell plates.

1- Seed cells of interest in a multi-well plate following standard protocols. (See an example of plate template in Figure 2)

2- Incubate the cells until 50-70% confluence is reached, approximately 10⁴ cells/well (times may vary depending on the cell type).

3- Prepare a 1:100 dilution of the non-fluorescent CONTROL SPAchips in cell culture medium. Mix thoroughly (vortex). Do not spin the tubes.

4-Dilute ASSAY SPAchip[®] stock solution in cell culture medium to obtain a final SPAchip[®] to cell ratio of 2:1. Alternatively, simply calculate 1 µL of stock per assay well and dilute 1:100 in cell culture medium. Vortex the ASSAY SPAchip[®] tube right before using it. Pipette up-and-down twice to homogenize the ASSAY SPAchip[®] solution and do not spin the tube.

! Optimization may be required depending on the cell type and assay conditions.

5- In the multi-well plate, aspire the cell medium and add 100 µL of CONTROL SPAchip[®] dilution to the control wells (see Figure 2). Vortex the tube right before use. Add 100 µL ASSAY SPAchip[®]-containing fresh culture medium per well. Homogenize the solution by pipetting up and down often.

6- Incubate overnight in a cell incubator to let the cells to internalize the SPAchip[®]. Internalization rate may vary depending on the cell subtype but should be over 25%.

⁴ For cells resistant to SPAchips internalization, lipofection may be necessary. Please, contact A4Cell for further information.

7- To include reference values, use some wells of the plate for calibrating the system (controls, calcium ionophores and/or calcium sequestrant such as Br-A23187 or BAPTA-AM in Figure 2). In such case, follow calibrator manufacturer's instructions.

8- Perform the experiment with your read-out platform. For long-time multiple-measurements assay (for example, monitorization during a week or month), keep the plate in proper conditions between each measurement and change the medium every 24-48 hours, depending on the cell subtype.




Figure 2: Example of 96-well plate configuration with positive controls and calibrators

c. SPAchip[®] Data Collection

- Image Acquisition:

1- Use transmitted light to localize at least three fields of interest per well (if using fluorescent light sources use long-wavelength light -561nm and over- in order to avoid photobleaching of the calcium-sensitive probes).

2- Use z-stacks to capture and set the best focus and the z-position of each SPAchip[®]. In confocal systems, reflected light can be used to visualize the SPAchips and the cells. Brightfield images or cell masking stains are recommended for intracellular quantification.

3- Adjust the acquisition parameters of the test probes: Ex 488/Em 520 (FITC fluorescence filter or similar) to avoid pixel saturation. If using calibrators, use A23187-treated wells to adjust this value.

4- Launch the experiment and save the images. Do not change settings parameters while acquisition.




Figure 4: a) Fluorescence signal of SPAchips in an imaging plate with cell culture medium. b) Representative image of HeLa cells with internalized SPAchips (red) stained with DAPI (blue) and phalloidin (green).

- Flow Cytometry:

- 1- Use FCS and SSC to locate SPAchips and cells independently by their size and complexity.
- 2- Create fluorescence histograms for green channel. Use cells treated with *CONTROL* SPAchips to set the fluorescence thresholds in both populations (free SPAchips and cells).
- 3- Analyze a tube of control cells treated with *ASSAY* SPAchips to verify the fluorescence threshold set in step 2 for gate “Cells”. Positive cells will appear as a separate peak in the histogram.

d. Quantitative analysis

1- **A)** Image Analysis: Select SPACHIP® fluorescent channel signal and measure the intensity of segmented Regions of Interest (ROIs) in every SPACHIP® with the image analysis software.

B) For flow cytometry assays: Analyze each sample with the adjusted gating parameters and quantify green channel detection (mean \pm SD or median \pm rSD).

2-Quantify the mean fluorescence of extra and intracellular SPACHIP® populations.

3- Export your data to a spreadsheet software.




Figure 5: Graph showing normalized relative fluorescence intensity values of CytoCHECK SPACHIP® Calcium Single Detection at different calcium concentration conditions.