

User Protocol

CytoCHECK SPAchip[®]

CALCIUM AND pH MULTI-DETECTION KIT

Product Reference: M-001-PC

For use with fluorescence microscopy, cell imaging platforms and flow cytometry.
Research use only.





CytoCHECK SPChip® CALCIUM and pH MULTI-DETECTION KIT

Product Reference: M-001-PC

For use with fluorescence microscopy, cell imaging platforms and flow cytometry.

1. General Introduction

CytoCHECK SPChip® assay kits are novel cell-based assays for living single-cell developed by A4cell that bring together the fields of nanotechnology and cell biology. CytoCHECK SPChip® kits are composed of fluorescently labeled silicon microparticles -SPChips®- that can be internalized in the cytosol of cultured cells to monitor changes in specific intracellular analyte concentrations for long periods of time.

CytoCHECK SPChip® Calcium and pH Multi-Detection Kit allows measuring cytosolic Calcium and intracellular and extracellular pH levels by changes in fluorescence intensity, which allows a more comprehensive study of the living single-cell physiology and maximizes the performance of most imaging analyzers. The main advantage of this product is the possibility of detecting Calcium and pH in the same SPChip.

CytoCHECK SPChip® Multi-Detection kits are optimized for its use with confocal microscopes and HCS/HCA analyzers with 20X or over magnification objectives, yet widefield fluorescence microscopes and imaging system with fixed wavelength filters can be also used. Flow cytometers have additionally been validated to analyze intra and extracellular SPChips.

SPChip® kits can be handled as common fluorophores and chemical probes for cell biology. After resuspending them in the assay buffer (a SPChip®-to-cell ratio of 2:1 is recommended), an overnight incubation is required to allow SPChip® to be incorporated in the cytosol¹. SPChip® will remain in the cytosol for longer than one month to monitor cell culture progression.

¹ CytoCHECK SPChip® kits are optimized for its use with adherent cells. Some cells may need a lipofection step to maximize SPChip internalization. Please contact with A4Cell for further information or questions about your cell model.



2.Contents and Storage

Each CytoCHECK SPChip® Calcium and pH Multi-Detection Kit contains:

Materials		Units
ASSAY SPChip® tube (embedded in a solid fluorescence-protective soluble film)	~2.5x10 ⁶ ASSAY SPChips	1
ASSAY buffer tube (sterile, cell culture suitable)	5 mL	1
CONTROL SPChip® tube (non-fluorescent, ready-to-use)	~2.5x10 ⁵ CONTROL SPChips/100 µL	1



Storage Information

Upon receipt, each kit should be stored protected from light at 2-8°C.

CytoCHECK SPChip® Calcium and pH Multi-Detection Kit Workflow

SPChip® Assay Workflow	Notes
<p style="writing-mode: vertical-rl; transform: rotate(180deg);">ASSAY SPChip® Preparation</p> <p style="text-align: center;">Dissolve ASSAY SPChip® tube film</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Wash with assay buffer } x2</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Resuspend</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Store</p>	<p>Thoroughly dissolve ASSAY SPChip® tube solid film in assay buffer (vortexing might be required).</p> <p>Centrifugate tubes at 4,300 xg for 10 min. Discard supernatant and keep SPChip® containing pellet.</p> <p>Resuspend the SPChips in 1mL ASSAY Buffer. Repeat x2.</p> <p>Resuspend the ASSAY SPChips in 100 µL ASSAY Buffer to obtain a stock solution of 2.5x10⁴ SPChips/µL*.</p> <p>ASSAY SPChip® stock solution can be stored protected from light at 4°C for up to 6 months.</p>
<p style="writing-mode: vertical-rl; transform: rotate(180deg);">Cell Assay</p> <p style="text-align: center;">Cell Culture</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">SPChip dilution</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Addition to Cell Culture</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Incubation</p>	<p>Seed cells (optimization will be required depending on the cell type and assay conditions) in a multi-well plate. Incubate the cells up to 50-70% confluence*.</p> <p>Dilute CONTROL and ASSAY SPChip stock solution in cell culture medium to obtain a final SPChip® to cell ratio of 2:1*.</p> <p>Aspire the cell medium and add SPChip®-loaded fresh culture medium*.</p> <p>Incubate overnight in a cell incubator (37°C, 5% CO₂) to allow SPChip® internalization by the cells.</p>
<p style="writing-mode: vertical-rl; transform: rotate(180deg);">Fluorescence Quantification</p> <p style="text-align: center;">Data Acquisition</p> <p style="text-align: center;">↓</p> <p style="text-align: center;">Data Analysis</p>	<p>Use green fluorescence channel to detect Calcium green signal (excitation 488 nm, emission 520 nm) and use red fluorescence channels to detect pH red ratiometric signal by exciting at 546 nm and collecting emission at 610 and 707 nm. Ratiometric intensities of both emission signals are required for optimal results.</p> <p>Analyze fluorescence intensity units and normalize the data to control conditions.</p>



*This cell culture protocol only provides a guideline, and it should be modified depending on specific needs. **Please, read whole protocol before starting.**

3. Materials to be Supplied by the User

- Pipettes and pipette tips (1-10 μ L, 2-20 μ L, 20-200 μ L, 200-1000 μ L).
- Multi-channel pipettes and pipette tips (50-300 μ L) (not mandatory but desirable).
- Vortex mixer.
- Mini-centrifuge.
- Reagent reservoirs.
- Cell culture conditions:
 - Cell culture facilities.
 - Cell culture plate (multi-well).
 - Cell culture media (phenol red free is highly recommended) according to specific cell line.
- Cell imaging system (HCS/HCS, Confocal microscope, imaging plate reader, etc.) with the appropriate fluorescence excitation and emission wavelength filters. At least 20x magnification objective is required for quantitative analysis.
- Image analysis software (Contact A4cell staff for support).

If quantitative analysis is required, it will be necessary:

- Intracellular pH calibrators.
- Calcium flux positive and negative controls (e.g. Br-A23187 and BAPTA-AM).

4. Assay Procedure

NOTE: Work under sterile conditions. Protect the tubes from light, specially from UV light.

a. ASSAY SPChip[®] Preparation:

1. Add 1 mL of assay buffer to ASSAY SPChip[®] tubes and mix until complete solubilization of the membrane (vigorous **vortexing** required).
2. Centrifuge the tubes at approx. 4,300 xg (8,000 RPM in a 6 cm minispin rotor) for 10 minutes.
3. Carefully, aspire and discard supernatant and wash pellets with 1 mL of solubilization medium (**Look out not to aspire the pellet!**). Repeat steps 1 and 2.
4. Carefully, aspire and discard supernatant. Resuspend the pellet in **100 µL of assay buffer** to obtain a SPChip[®] stock solution. This results in approximately **2.5x10⁴ SPChips/µl¹**. Once prepared, ASSAY SPChip[®] tube can be stored at 2-8°C protected from light for up to 6 months.



It is essential to properly dissolve the SPChip[®] solid film for the success of the assay.

b. ASSAY PROTOCOL (standard protocol for 96-well plate to work with adherent cells):

NOTE: ASSAY and CONTROL SPChip[®] dilutions suggested in this section are suggested for 96-well plates. Optimization might be required for different cell types or multiwell plates.

Goal: To demonstrate intracellular pH and Calcium monitorization in adherent cells under different treatment conditions by using SPChip[®] nanodevices with fluorescence microscopy techniques.

1. **Seed cells** of interest in a multi-well plate following standard protocols. See an example of plate template in **Figure 1**: control wells are non-fluorescent SPChips; *pH 4.5, 5.5, 6.5 and 7.5* are quality control wells for calculating an intracellular pH calibration curve where measurements should be interpolated; Br-A23187 and BAPTA-AM are two examples of Calcium ionophore and Calcium chelator that can be used as positive controls for increasing and decreasing Calcium signal respectively. Blank conditions with the sample medium should be used as negative control; *sample 1, 2, etc* are the measurements of interest whose values should be interpreted according to the Calcium controls and pH calibrators.

2. Incubate the cells until 50-70% confluence is reached, approximately 10⁴ cells/well (times may vary depending on the cell type)².

3. Prepare a 1:100 dilution of the non-fluorescent CONTROL SPChips in cell culture medium³. Mix thoroughly (vortex). Do not spin the tubes.



¹It is highly recommended the previous counting of SPChips in a Neubauer chamber. Check SPChip[®] count using a routine brightfield microscope.



²Optimization may be required depending on the cell type and assay conditions. For cells resistant to SPChips internalization, lipofection may be necessary. Please, contact A4cell for further information.



³We recommend to previously test SPChip[®] internalization by performing a preliminary experiment using CONTROL SPChips and checking for % of internalization using a brightfield microscope. This ensures the appropriate quantity and internalization of SPChips, allowing cell assay optimization.

4. **Dilute** ASSAY SPACHIP® stock solution in cell culture medium to obtain a final SPACHIP® to cell ratio⁴ of 2:1. Alternatively, simply calculate 1 µL of stock per assay well and dilute 1:100 in cell culture medium. Vortex the ASSAY SPACHIP® tube right before using it. Pipette up-and-down twice to homogenize the ASSAY SPACHIP® solution and do not spin the tube.
5. If fluorescence baseline is required, CONTROL SPACHIP® samples (non-fluorescent) can be included in assay design. In a 96-well plate, aspire the cell medium and add 100 µL of CONTROL SPACHIP® **dilution** to the control wells (Figure 1). Vortex the tube right before use. Add 100 µL ASSAY SPACHIP®-containing fresh culture medium per well⁵. Homogenize the solution by pipetting up and down several times.
6. **Incubate** overnight in a cell incubator to let the cells to internalize the SPACHIP®. Internalization rate may vary for each cell type, regular values for adherent cells range from 25 to 50% depending on confluence, SPACHIP®-to-cell ratio and assay duration.
7. If quantification is needed, use some wells of the plate for calibrating the system (controls and intracellular pH calibrators in the columns 1-4, see example of plate template below in Figure 1). In such case, follow calibrator manufacturer's instructions.
8. Perform the experiment with your read-out platform. For long-time kinetic assay (for example, monitorization during one week or month), keep the plate under proper conditions (CO₂, T, humidity) between each measurement and change the medium accordingly, depending on the cell type.




Figure 1: Example of 96-well plate configuration with positive controls and calibrators.



⁴It is highly recommended to test SPACHIP® internalization in the specific cell line. The appropriate SPACHIP® to cell ratio might vary depending on the cell type and assay conditions. Thus, we recommend to previously optimize SPACHIP® to cell ratio in previous step, using CONTROL SPACHIPs.



⁵Phenol Red can quench the signal of some visible wavelength fluorescence dyes and cause problems when used during live cell imaging, interfering in the read-out. We recommend to add a fresh media that does not contain Phenol Red in the formulation.

c. SPACHIP® Data Collection

Image Acquisition

1. Set objective to 20x magnification or over (enough resolution for 3 μm SPACHIP®)
2. Use transmitted light to localize the region of interest (if using fluorescent light sources use long-wavelength light -561 nm and over- in order to avoid photobleaching of the Calcium and pH sensitive probes).
3. Use Z-stacks to capture and set the best focus and the Z position of each SPACHIP®. Brightfield images or cell masking stains are recommended for intracellular quantification.
4. Adjust the acquisition parameters for SPACHIP® fluorescent signals. Calcium channel corresponding to Ex 488 nm /Em 520 nm (FITC fluorescence filter or similar) and pH: Ex 546 nm /Em 610 & 707 nm. If using calibrators or positive controls, to avoid pixel saturation, use higher fluorescence intensity conditions corresponding to Br-A23187 for Calcium and pH 7.5 wells to adjust this value.
5. Launch the experiment and save the images. Do not change settings parameters while acquisition.




Figure 2: **A)** Fluorescence signal of SPACHIPs in an imaging plate with cell culture medium. Representative images of SH-SY5Y cells with internalized multi-detection SPACHIPs for **B)** Calcium signal exciting at 488 nm and emitting at 520 nm (green); and **C)** pH signal exciting at 546 nm and collecting emission at 610 and 707 nm (red). Images B) and C) correspond to the same field but different fluorescence channels. Nuclei are stained in blue for images B) and C).

Flow Cytometry

1. Use FCS and SSC to locate SPACHIPs and cells independently by their size and complexity.
2. Create fluorescence histograms for green and red channels. Use cells treated with CONTROL SPACHIPs to set the fluorescence thresholds in both populations (free SPACHIPs and cells).
3. Analyze a tube of control cells treated with ASSAY SPACHIPs to verify the fluorescence threshold set in step 2 for gate “Cells”. Positive cells will appear as a separate peak in the histogram.

d. Quantitative Analysis

1. **Image Analysis:** select SPChip® fluorescent channel signals (emission at 610 and 707 nm for pH red signal and emission at 520 nm for calcium green signal) and measure the intensity of segmented Regions of Interest (ROIs) in every SPChip® with the image analysis software. Quantify the mean fluorescence of extra and intracellular SPChip® populations. Afterwards, ratio of these pH measurements (emission at 610 and 707 nm) can be taken to increase signal to noise ratio.

- For flow cytometry assays: Analyze each sample with the adjusted gating parameters and quantify two channel detections (mean \pm SD or median \pm rSD).

2. **Quantify** the mean fluorescence of extra and intracellular SPChip® populations.

3. **Export** your data to a spreadsheet software.

4. **To quantify** intracellular **pH**, use the fluorescence of the calibration wells to plot ratiometric fluorescence units vs calibration pH. Interpolate the values of the sample wells which fits to a sigmoidal curve (Figure 3A) to obtain the intracellular pH value of the sample wells. It is important to interpolate sample values in the intracellular calibration curves instead of extracellular media.

5. **To quantify** extracellular **calcium** levels, the fluorescence of buffers with different calcium concentrations can be measured, obtaining a calibration curve (Figure 3B). To determine the variation of intracellular Calcium levels, use ionophores and/or Calcium chelators to obtain relative values to reference control conditions.




Figure 3: CytoCHECK SPChip® Calcium and pH Multi-Detection Kit calibration curve. **a)** Graph showing ratiometric normalized mean fluorescence intensity (MFI) values at different pH. Ratiometric values were obtained by dividing $\lambda_{em2}/\lambda_{em1}$ emission signals in HCS-Operetta equipment with the excitation in the range $\lambda_{exc}=545/15$ nm. **b)** Extracellular Calcium calibration curve.

* Contact A4cell staff for assistance with imaging setup and data analysis. Ask for our **cell biology image analysis guide** at labservices@a4cell.com .